

Product Datasheet

Mouse Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit (orb1668191)

Description	Mouse Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-4000 ng/mL. Sensitivity: 19.3 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.