

Product Datasheet

General 8-Hydroxydeoxyguanosine (8-OHdG) ELISA Kit (orb1668192)

Description	General 8-Hydroxydeoxyguanosine (8-OHdG) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-4000 pg/mL. Sensitivity: 18.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.