

Product Datasheet

Classical Swine Fever Virus Anitibody ELISA Kit (orb1668210)

Description	Classical Swine Fever Virus Anitibody ELISA Kit
Reactivity	Virus
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.