

Product Datasheet

Human Hepatitis B virus pre S1 antigen (HBV preS1-Ag) ELISA kit (orb1668232)

Description	Human Hepatitis B virus pre S1 antigen (HBV preS1-Ag) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.