

Product Datasheet

Pig Pseudorabies virus IgG antibody (PRV-Ab-IgG) ELISA Kit (orb1668332)

Description	Pig Pseudorabies virus IgG antibody (PRV-Ab-IgG) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.