

Product Datasheet

Beclin 1 (phospho-S64) antibody (orb166835)

Description	Rabbit polyclonal antibody to Beclin 1 (phospho-S64)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	BECN1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51896
Uniprot ID	Q14457
NCBI	NP_003757.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

