

Product Datasheet

CCNB2 (phospho-S22) antibody (orb166837)

Description Rabbit polyclonal antibody to CCNB2 (phospho-S22).

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications DOT

Immunogen Synthetic Peptide

Target CCNB2

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 45282

Uniprot ID [O95067](#)

NCBI [NP_004692.1](#)

Dilution Range DOT: 1:500

Expiration Date 12 months from date of receipt.

