

Product Datasheet

Human Lipoxin A4 (LXA4) ELISA Kit (orb1668411)

Description	Human Lipoxin A4 (LXA4) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 493.8-40000 pg/mL. Sensitivity: 175.4 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.