

Product Datasheet

Human Trypsinogen-2 (TRY2) ELISA Kit (orb1668419)

Description Human Trypsinogen-2 (TRY2) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.781-50 ng/mL.
Sensitivity: 0.28 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.