

Product Datasheet

Human Soluble fibrin monomer complex (SFMC) ELISA Kit (orb1668465)

Description	Human Soluble fibrin monomer complex (SFMC) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.125-200 µg/mL. Sensitivity: 1 µg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.