

Product Datasheet

ULK1 (phospho-S467) antibody (orb166857)

Description	Rabbit polyclonal antibody to ULK1 (phospho-S467).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	ULK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	112631
Uniprot ID	O75385
NCBI	NP_003556.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

