

Product Datasheet

Rat Low Density Lipoprotein cholesterol (HDL-C) ELISA Kit (orb1668608)

Description	Rat Low Density Lipoprotein cholesterol (HDL-C) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 ng/mL. Sensitivity: 2.71 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.