

Product Datasheet

CCNB3 (phospho-S1063) antibody (orb166861)

Description	Rabbit polyclonal antibody to CCNB3 (phospho-S1063).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	Ccnb3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	158969
Uniprot ID	Q810T2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

