

Product Datasheet

p27Kip1 (phospho-S83) antibody (orb166863)

Description	Rabbit polyclonal antibody to p27Kip1 (phospho-S83).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	CDKN1B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22073
Uniprot ID	P46527
NCBI	NP_004055.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

