

Product Datasheet

Monkey Fibrinogen degradation product (FDP) ELISA Kit (orb1668637)

Description	Monkey Fibrinogen degradation product (FDP) ELISA Kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.