

Product Datasheet

Hamster Follicle stimulating hormone (FSH) ELISA Kit (orb1668643)

Description	Hamster Follicle stimulating hormone (FSH) ELISA Kit
Reactivity	Hamster
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.