

Product Datasheet

Bovine Fibroblast growth factor 21 (FGF21) ELISA Kit (orb1668662)

Description	Bovine Fibroblast growth factor 21 (FGF21) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.