

Product Datasheet

LRP2 (Megalin) antibody (orb166882)

Description	Rabbit polyclonal antibody to LRP2 (C-term) (Megalin).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	LRP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	521958
Uniprot ID	P98164
NCBI	NP_004516.2
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.

