



Product Datasheet

NRG2 antibody (orb166883)

Description	Rabbit polyclonal antibody to NRG2 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	NRG2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	91679
Uniprot ID	O14511
NCBI	NP_001171864.1 , NP_053586.1 , NP_053584.1 , NP_004874.1 , NP_053585.1
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



WB analysis
of HT-29
whole cell
lysate u...