

Product Datasheet

Human Anti-rubella virus IgG antibody (RV-Ab-IgG) ELISA Kit (orb1668862)

Description	Human Anti-rubella virus IgG antibody (RV-Ab-IgG) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Detection antibody. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.