

Product Datasheet

Human Anti-mutated citrullinated vimentin antibody (MCV-Ab) ELISA Kit (orb1668868)

Description Human Anti-mutated citrullinated vimentin antibody (MCV-Ab) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 5 U/mL-400 U/mL. Sensitivity: 3.25 μ /mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.