

Product Datasheet

Horse Brain-derived neurotrophic factor (BDNF) ELISA Kit (orb1668927)

Description	Horse Brain-derived neurotrophic factor (BDNF) ELISA Kit
Reactivity	Equine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.25-2000 pg/mL. Sensitivity: 12.6 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.