

## Product Datasheet

### Human Muscarinic acetylcholine receptor M3 (CHRM3) ELISA Kit (orb1668959)

|                            |                                                                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Muscarinic acetylcholine receptor M3 (CHRM3) ELISA Kit                                                                                                                           |
| <b>Reactivity</b>          | Human                                                                                                                                                                                  |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                                                  |
| <b>Application notes</b>   | Detect Range: 0.78-50 ng/mL. Sensitivity: 0.30 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                        |