

Product Datasheet

Myoglobin antibody (orb166899)

Description	Mouse monoclonal antibody to Myoglobin.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MB
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	This antibody is purified through a protein G column, followed by dialysis against PBS.
MW	17184 Da
Uniprot ID	P02144
NCBI	NP_976312.1, NP_005359.1, NP_976311.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

