

Product Datasheet

Erk2 antibody (orb166903)

Description	Mouse monoclonal antibody to Erk2.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MAPK1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	41390
Uniprot ID	P28482
NCBI	NP_620407.1 , NP_002736.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

