

Product Datasheet

Mouse Nuclear factor erythroid 2-related factor 3 (NFE2L3) ELISA Kit (orb1669120)

Description	Mouse Nuclear factor erythroid 2-related factor 3 (NFE2L3) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.057 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.