

Product Datasheet

p53 antibody (orb166919)

Description	Mouse monoclonal antibody to p53
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	TP53
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	43653
Uniprot ID	P04637
NCBI	NP_001119586.1, NP_001119590.1, NP_001119585.1, NP_001263625.1, NP_001263626.1, NP_001119584.1, NP_001263624.1, NP_00126, NP_001119589.1, NP_001263627.1, NP_001119587.1, NP_000537.3, NP_001263628.1, NP_001119588.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

