

Product Datasheet

Histone H3 (phospho-S10) antibody (orb166933)

Description	Mouse monoclonal antibody to Histone H3 (phospho-S10).
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	H3-3A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	15328
Uniprot ID	P84243
NCBI	NP_005315.1, NP_002098.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

