

Product Datasheet

Mouse Dipeptidyl peptidase 9 (DPP9) ELISA Kit (orb1669357)

Description	Mouse Dipeptidyl peptidase 9 (DPP9) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 25 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.