

Product Datasheet

PRDM10 antibody (orb166939)

Description	Mouse monoclonal antibody to PRDM10.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	PRDM10
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	130136
Uniprot ID	Q9NQV6
NCBI	NP_955469.1 , NP_955470.1 , NP_064613.2 , NP_955471.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

