

Product Datasheet

Fish Interleukin 1 β (IL-1 β) ELISA Kit (orb1669398)

Description	Fish Interleukin 1 β (IL-1 β) ELISA Kit
Reactivity	Fish
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.125-200 pg/mL. Sensitivity: 1.1 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.