

Product Datasheet

Human Inositol 1,4,5,-trisphosphate (IP3) ELISA Kit (orb1669400)

Description	Human Inositol 1,4,5,-trisphosphate (IP3) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 123.5-10000 pg/mL. Sensitivity: 49.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.