

Product Datasheet

SUMO1 antibody (orb166942)

Description	Mouse monoclonal antibody to SUMO1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	SUMO1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	11557
Uniprot ID	P63165
NCBI	NP_001005782.1 , NP_001005781.1 , NP_003343.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

