



Product Datasheet

UBC9 antibody (orb166944)

Description	Mouse monoclonal antibody to UBC9.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	UBE2I
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	18007
Uniprot ID	P63279
NCBI	NP_919236.1, NP_919237.1, NP_003336.1, NP_919235.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse spleen tissue.