

Product Datasheet

Rat Prostaglandin E2 (PG-E2) ELISA Kit (orb1669453)

Description Rat Prostaglandin E2 (PG-E2) ELISA Kit

Reactivity Rat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 3.12-200 pg/mL.
Sensitivity: 1.04 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.