




Product Datasheet

NUCL antibody (orb166948)

Description	Mouse monoclonal antibody to NUCL.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	NCL
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1?
Clonality	Monoclonal
MW	76614
Uniprot ID	P19338
NCBI	NP_005372.2
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...