

Product Datasheet

Human Ubiquitin carboxyl-terminal hydrolase BAP1 (BAP1) ELISA Kit (orb1669502)

Description Human Ubiquitin carboxyl-terminal hydrolase BAP1 (BAP1) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 78.1-5000 pg/mL. Sensitivity: 30 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.