

Product Datasheet

Hamster Very low density lipoprotein (VLDL) ELISA Kit (orb1669508)

Description	Hamster Very low density lipoprotein (VLDL) ELISA Kit
Reactivity	Hamster
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8-500 ng/mL. Sensitivity: 3.8 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.