




Product Datasheet

NGFR antibody (orb166956)

Description	Mouse monoclonal antibody to NGFR
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	NGFR
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2b,k
Clonality	Monoclonal
MW	45183
Uniprot ID	P08138
NCBI	NP_002498.1
Dilution Range	IF/ICC: 1:100, WB: 1:500-1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart tis...



Immunohistochemical analysis of formalin...