

Product Datasheet

Human Collagen Type III (Col III) ELISA Kit (orb1669571)

Description	Human Collagen Type III (Col III) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.5-800 ng/mL. Sensitivity: 6.4 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.