



Product Datasheet

LAMP2 antibody (orb166960)

Description	Mouse monoclonal antibody to LAMP2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant Protein
Target	LAMP2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	44961
Uniprot ID	P13473
NCBI	NP_001116078.1 , NP_054701.1 , NP_002285.1
Dilution Range	WB: 1:100-1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...