




Product Datasheet

RPS6KB2 antibody (orb166961)

Description	Mouse monoclonal antibody to RPS6KB2.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	RPS6KB2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	53455
Uniprot ID	Q9UBS0
NCBI	NP_003943.2
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse NIH-3T3 c...



WB analysis of Jurkat MCF-7 mouse NIH/3T...