

Product Datasheet

THEM2 antibody (orb166966)

Description	Mouse monoclonal antibody to THEM2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	ACOT13
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	14960
Uniprot ID	Q9NPJ3
NCBI	NP_001153566.1 , NP_060943.1
Dilution Range	WB: 1:100-1000
Expiration Date	12 months from date of receipt.

