

Product Datasheet

LXN antibody (orb166976)

Description	Mouse monoclonal antibody to LXN.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	LXN
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	25750
Uniprot ID	Q9BS40
NCBI	NP_064554.3
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

