



Product Datasheet

PRKAA2 antibody (orb166983)

Description	Mouse monoclonal antibody to PRKAA2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant Protein
Target	PRKAA2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2a,K
Clonality	Monoclonal
MW	62320
Uniprot ID	P54646
NCBI	NP_006243.2
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of CEM cell line lysate.