

## Product Datasheet

### Human Beta-defensin 130 (DEFB130) ELISA Kit (orb1669856)

|                            |                                                                                                                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Beta-defensin 130<br>(DEFB130) ELISA Kit                                                                                                                                                                          |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                   |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                   |
| <b>Note</b>                | For research use only                                                                                                                                                                                                   |
| <b>Application notes</b>   | Detect Range: Request Information.<br>Sensitivity: Request Information.<br>Assay Time: 1-4.5 h. Detection<br>wavelength: 450 nm. Sample<br>Volume: 1-200 µL. Sample Type:<br>Serum, Plasma, Other biological<br>fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                         |