

Product Datasheet

EIF4E2 antibody (orb166986)

Description	Mouse monoclonal antibody to EIF4E2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	EIF4E2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	28362
Uniprot ID	O60573
NCBI	NP_004837.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

