

Product Datasheet

ATG4B antibody (orb166987)

Description	Mouse monoclonal antibody to ATG4B.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	ATG4B
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	44294
Uniprot ID	Q9Y4P1
NCBI	NP_847896.1 , NP_037457.3
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

