



Product Datasheet

PTPN11 antibody (orb166998)

Description	Mouse monoclonal antibody to PTPN11.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	PTPN11
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	68011
Uniprot ID	Q06124
NCBI	NP_002825.3 , NP_542168.1
Dilution Range	WB: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...