



Product Datasheet

CSNK2B antibody (orb166999)

Description	Mouse monoclonal antibody to CSNK2B.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	CSNK2B
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	24942
Uniprot ID	P67870
NCBI	NP_001269314.1, NP_001311.3
Dilution Range	WB: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse NIH-3T3 c...